


SyndoraAlto.



POWERED BY




altoEnhance

Integrated with **Silver Shield**

Basic Product Information

v4.2





altoEnhance is an intelligent system which sends the alerts from your Nurse Call System to the android devices carried by care staff.

The Nurse Call Alerts App allows care staff to accept and efficiently respond to alerts made from residents, by showing the alerts to staff who are not in visible range of the call bell screen panels, and ultimately reduce the noise within the home.



Leading, AI-powered Fall Detection & Prevention






- Eliminate false alerts with AI sensing**
- No cameras, mics, or wearables**
- Plug-and-play, w/o learning period**
- Over 97% accuracy**
- Differentiates between people and animals**
- 100% accurate room presence and real-time location view**
- Unparalleled 6x6m coverage, no dead zones**
- Tracks habits for sharper insights and helps prevent bed exit falls**
- Wellbeing Monitoring incl. breathing rate**


Pontosense Silver Shield is an intelligent monitoring solution that uses advanced radar technology to detect falls, breathing patterns, and presence within a resident's room.

The system provides real-time insights and silent alerts to care staff, helping them respond quickly to incidents without the need for wearable devices or intrusive cameras.

Silver Shield enhances safety, supports proactive care, and reduces disturbance within the home by delivering accurate, discreet monitoring around the clock.




How does Silver Shield and altoEnhance work together?



Falls, Out of Bed and Room Exit events are seamlessly integrated to **altoEnhance**, which appear in a familiar fashion as regular nurse call alerts.

This allows staff to accept and respond to **Silver Shield** events directly within the nurse call alerts app which they are already using - no additional training required.

How are Falls from Silver Shield delivered to the Nurse Call Alerts App?



Falls, Out of Bed and Room Exit events are sent directly to the carers devices.

A staff member accepts the alert, and attends to the location to reset the alert and provide any necessary care.



01635 019608



hello@spark-care.co.uk



spark-care.co.uk